



96-well microplates BRANDplates® with C-bottom

Polystyrene, untreated or optimized for immunological and cell culture applications.

Cat.No.	Bottom	Surface	Color	Type	Sterile	max. Vol.	Lid
7816 03	C	pureGrade™	transparent	standard	-	350 µl	-
7816 63	C	pureGrade™ S	transparent	standard	+	350 µl	+
7817 23	C	immunoGrade™	transparent	standard	-	350 µl	-
7817 26	C	immunoGrade™	white	standard	-	350 µl	-
7817 29	C	immunoGrade™	black	standard	-	350 µl	-
7817 83	C	hydroGrade™	transparent	standard	-	350 µl	-
7817 86	C	hydroGrade™	white	standard	-	350 µl	-
7817 89	C	hydroGrade™	black	standard	-	350 µl	-
7818 43	C	lipoGrade™	transparent	standard	-	350 µl	-
7818 46	C	lipoGrade™	white	standard	-	350 µl	-
7818 49	C	lipoGrade™	black	standard	-	350 µl	-
7819 63	C	cellGrade™	transparent	standard	+	350 µl	+
7819 06	C	inertGrade™	white	standard	+	350 µl	+
7819 09	C	inertGrade™	black	standard	+	350 µl	+

untreated
immunological
cell culture

All BRANDplates® microplates are manufactured according to the requirements of the ANSI-/SLAS-standards 1 to 4 - 2004. This ensures the compatibility of the BRANDplates® when performing automated processes. Subject to technical modification without notice. Errors excepted.