



## 96-well microplates BRANDplates<sup>®</sup> with U-bottom

Polystyrene, untreated or optimized for immunological and cell culture applications.

Cat.No.	Bottom	Surface	Color	Type	Sterile	max. Vol.	Lid
7816 00	U	pureGrade™	transparent	standard	-	330 µl	-
7816 04	U	pureGrade™	white	standard	-	330 µl	-
7816 07	U	pureGrade™	black	standard	-	330 µl	-
7816 60	U	pureGrade™ S	transparent	standard	+	330 µl	+
7817 20	U	immunoGrade™	transparent	standard	-	330 µl	-
7817 24	U	immunoGrade™	white	standard	-	330 µl	-
7817 27	U	immunoGrade™	black	standard	-	330 µl	-
7817 80	U	hydroGrade™	transparent	standard	-	330 µl	-
7817 84	U	hydroGrade™	white	standard	-	330 µl	-
7817 87	U	hydroGrade™	black	standard	-	330 µl	-
7818 40	U	lipoGrade™	transparent	standard	-	330 µl	-
7818 44	U	lipoGrade™	white	standard	-	330 µl	-
7818 47	U	lipoGrade™	black	standard	-	330 µl	-
7819 60	U	cellGrade™	transparent	standard	+	330 µl	+
7819 00	U	inertGrade™	transparent	standard	+	330 µl	+
7819 04	U	inertGrade™	white	standard	+	330 µl	+
7819 07	U	inertGrade™	black	standard	+	330 µl	+

untreated
immunological
cell culture

All BRANDplates<sup>®</sup> microplates are manufactured according to the requirements of the ANSI-/SLAS-standards 1 to 4 - 2004. This ensures the compatibility of the BRANDplates<sup>®</sup> when performing automated processes. Subject to technical modification without notice. Errors excepted.