



96-well strip plates BRANDplates[®] with F-bottom

Polystyrene, untreated or optimized for immunological and cell culture applications.

Cat.No.	Bottom	Surface	Color	Type	Sterile	max. Vol.	use of single wells
7823 00	F	pureGrade [™]	transparent	strip plate	-	360 µl	-
7823 01	F	pureGrade [™]	transparent	strip plate*	-	360 µl	+
7823 05	F	immunoGrade [™]	transparent	strip plate	-	360 µl	-
7823 06	F	immunoGrade [™]	transparent	strip plate*	-	360 µl	+

* 8-well Strips in single well holding frame.

untreated
immunological
cell culture

All BRANDplates[®] microplates are manufactured according to the requirements of the ANSI/SLAS-standards 1 to 4 - 2004. This ensures the compatibility of the BRANDplates[®] when performing automated processes. Subject to technical modification without notice. Errors excepted.

