



desde 1967
AVÁNTIKA

MEZCLADORES VORTEX VELP

TABLA DE COMPARACIÓN



RX3



ZX3



ZX4

Stirring speed rpm		3000, constant	0 3000	0 3000
Orbital diameter mm		4.5	4.5	4.5
Operating mode	touch	•	•	
	continuous		•	•
	IR sensor			•
Digital display and timer				
Dimensions (WxHxD) mm		150x130x165	150x130x165	150x130x165
	(in)	(5.9x5.1x6.5)	(5.9x5.1x6.5)	(5.9x5.1x6.5)
Weight	Kg	2.7	2.7	2.7
	(lb)	(5.9)	(5.9)	(5.9)
Power Supply		100 240 V	100 240 V	100 240 V



TX4



CLASSIC



WIZARD

Stirring speed rpm		0 3000	0 3000	0 3000
Orbital diameter mm		4.5	4.5	4.5
Operating mode	touch		•	
	continuous	•	•	•
	IR sensor	•		•
Digital display and timer				
Dimensions (WxHxD) mm		150x130x165	180x70x220	180x70x220
	(in)	(5.9x5.1x6.5)	(7.1x2.8x8.7)	(7.1x2.8x8.7)
Weight	Kg	2.7	2.2	2.2
	(lb)	(5.9)	(4.9)	(4.9)
Power Supply				