



BRAND



Compatibility of Pipette Tips and Filter Tips from BRAND

	Pipette tips (µl)						Filter tips (µl)						
	0.1-20	0.5-20	1-50	2-200	5-300	50-1000	0.1-1	0.5-10	1-20	2-20	5-100	5-200	50-1000
BRAND®													
Transferpette® S													
0.1 - 1 µl	●						●						
0.1 - 2.5 µl	●	●	●					●	●	●			
0.5 - 10 µl	●	●	●				●	●	●				
2 - 20 µl				●	●					●	●		
10 - 100 µl				●	●					●	●		
20 - 200 µl				●	●					●	●	●	
100 - 1000 µl						●							●
Transferpette® electronic													
0.5 - 10 µl	●	●	●				●	●	●				
2 - 20 µl	●	●	●				●	●	●				
20 - 200 µl				●	●					●	●	●	
100 - 1000 µl						●							●
Eppendorf®													
Reference®													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●	●					●	●				
2 - 20 µl (grey)	●	●							●				
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●						●	●	●	
50 - 200 µl				●		●					●	●	
100 - 1000 µl						●							●
Reference® 2													
0.1 - 2.5 µl	●	●	●				●	●	●				
0.5 - 10 µl	●	●	●				●	●	●				
2 - 20 µl (grey)	●	●	●					●	●				
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●	●					●	●	●	
20 - 200 µl				●	●					●	●	●	
30 - 300 µl				●	●						●	●	
100 - 1000 µl						●							●
Research®													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●						●					
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●						●	●	●	
20 - 200 µl				●							●	●	
100 - 1000 µl						●							●
Research® plus													
0.1 - 2.5 µl	●	●						●					
0.5 - 10 µl	●	●	●					●	●				
2 - 20 µl (grey)	●	●	●						●				
2 - 20 µl (yellow)				●						●	●	●	
10 - 100 µl				●						●	●	●	
20 - 200 µl				●							●	●	
30 - 300 µl					●							●	
100 - 1000 µl						●							●
Xplorer®													
0.5 - 10 µl	●	●	●					●	●				
5 - 100 µl				●							●		
15 - 300 µl												●	
50 - 1000 µl						●							
Xplorer® plus													
0.5 - 10 µl	●	●	●					●	●				
5 - 100 µl				●							●		
15 - 300 µl												●	
50 - 1000 µl						●							●
BIOHIT®													
mLINE													
0.1 - 3 µl	●	●	●					●	●				
0.5 - 10 µl	●	●						●					
2 - 20 µl				●					●	●			
10 - 100 µl				●						●	●		
20 - 200 µl				●									
100 - 1000 µl						●							●

● = compatible □ = no information

Continued on next page

Compatibility of Pipette Tips and Filter Tips from BRAND

	Pipette tips (µl)					Filter tips (µl)							
	0.1-20	0.5-20	1-50	2-200	5-300	50-1000	0.1-1	0.5-10	1-20	2-20	5-100	5-200	50-1000
GILSON®													
Pipetman P													
0.2 - 2 µl	●	●	●					●	●				
1 - 10 µl	●	●	●					●	●				
2 - 20 µl				●						●			
20 - 100 µl				●							●	●	
50 - 200 µl				●								●	
200 - 1000 µl						●							●
Pipetman Neo													
0.2 - 2 µl	●	●	●					●					
1 - 10 µl	●	●	●										
2 - 20 µl				●						●			
10 - 100 µl				●							●	●	
20 - 200 µl				●								●	
100 - 1000 µl						●							●
Thermo													
FINNPIPETTE® Digital													
0.2 - 2 µl (micro)													
0.5 - 10 µl (micro)													
0.5 - 10 µl	●	●						●	●				
2 - 20 µl													
5 - 50 µl				●							●	●	
10 - 100 µl													
20 - 200 µl				●								●	
100 - 1000 µl						●							●
FINNPIPETTE® Focus													
0.3 - 3 µl	●	●	●										
0.5 - 5 µl													
1 - 10 µl (pink)	●	●						●			●		
1 - 10 µl (yellow)				●						●			
FINNPIPETTE® Focus													
3 - 30 µl (turquoise)													
3 - 30 µl (yellow)													
5 - 50 µl (turquoise)													
5 - 50 µl (yellow)													
10 - 100 µl				●							●	●	
30 - 300 µl					●								●
100 - 1000 µl						●							●
Abimed													
Discovery													
0.5 - 10 µl	●	●						●	●				
2 - 20 µl				●						●			
10 - 100 µl													
20 - 200 µl				●								●	
100 - 1000 µl						●							●
Labmate													
0.5 - 10 µl	●	●	●					●	●				
2 - 20 µl				●						●			
10 - 100 µl				●							●		
20 - 200 µl				●								●	
100 - 1000 µl						●							●

● = compatible □ = no information

Status as of: 1014

This table is based on general experiences and therefore we cannot assume any warranty. The customer should always check the compatibility under his working conditions.

Transferpette®, BRAND® as well as the logo design mark depicted here are trademarks of BRAND GMBH + CO KG, Germany. The following brands and trademarks as referenced in the current table are property of their respective third-party owners as listed below:

Referenced brand/trademark	Owner
BIOHIT®	Sartorius Biohit Liquid Handling Oy, Finland
Eppendorf®	Eppendorf AG, Germany
FINNPIPETTE®	Thermo Fisher Scientific Oy, Finland
GILSON®	Gilson, Inc., USA
Reference®	Eppendorf AG, Germany
Research®	Eppendorf AG, Germany
Xplorer®	Eppendorf AG, Germany